

PowerTek™

CAPITAL SYSTEM

Featuring...

PowerTek™ 
HIGH RES 3CCD DIGITAL VIDEO CAMERA

PowerTek™ II Plus
HIGH SPEED SHAVER SYSTEM

PowerPump® 
SYSTEM

BIOMET®
SPORTS MEDICINE





PowerTek™ HR

HIGH RES 3CCD DIGITAL VIDEO CAMERA

- 3-chip all digital camera with programmable buttons on compact camera head
- Unique software enhances depth perception for a natural color reproduction
- 300 watt Xenon light source with specially designed lenses for optimal light distribution and high color temperature
- Medical grade Sony 21" LCD monitor
- Sony medical grade DVD recorder offers up to three hours of recording time
- Sony color printer provides a superior print quality and color reproduction with resolution of approximately 400 dots per inch



**3CCD Camera
with Programmable Buttons**



**Medical-Grade Sony
21" LCD Monitor**



Sony Color Printer



Sony Medical DVD Recorder

For complete instructions, risks and warnings, see the appropriate operator's manual for each piece of equipment.



PowerTek™ II Plus

HIGH SPEED SHAVER SYSTEM

- Capable of speeds up to 12,000 RPM for efficient bone removal and meniscal cutting
Note: Blade use is not recommended at speeds above 5,000 RPM.
- High rotational force handpiece/high torque
- Multifunctional buttons on handpiece allows user to control functions from sterile field
- Full line of burs and blades available
- Quick release latch to rapidly change from blade to bur
- Tool mode allows rapid increase and decrease between speeds
- Ergonomic suction valve allows surgeons linear suction control for debris aspiration



High Speed Suction Control Handpiece

PowerPump®

SYSTEM

- Provides constant visualization by providing constant dilation, even when there are pressure changes in the joint
- Double measuring device incorporates two pressure sensors which may increase patient safety
- Offers an automatic software-controlled instrument recognition system which communicates with different brands of endoscopes by measuring the flow resistance to the instrument
- Measures flow and pressure
- Quick and easy to set up—just set the pressure and begin!



Ordering Information

Sony® 21" Flat Screen Monitor

906083

Sony® DVO Video Capture

906084

Sony® DVO Video Capture Footswitch

906089

Sony® Color Printer

906085

PowerTek™ Camera Console with Handpiece

906070

PowerTek™ Camera Handpiece

906072

Xenon Light Source

906068

Universal Light Cable

906117

Arthroscope

900203 4mm 30°
900204 4mm 70°
900123 2.7mm 30°

Endovideoscope

900205 4mm 30°
900206 4mm 70°

5.5mm Cannula (for use with 4mm Scope)

900006 Two Stopcock
900007 One Stopcock

4.1mm Cannula (for use with 2.7mm Scope)

900026 Two Stopcock
900027 One Stopcock

5.5mm Quick-Release Cannula (for use with 4mm Scope)

900057 One Stopcock

Trocar

900011 5.5mm
900059 5.5mm Quick Release
900034 4.1mm

Obturator

900016 5.5mm
900058 5.5mm Quick-Release
900035 4.1mm

Coupler

907412 18mm
907413 22mm
907416 24mm
907414 25mm

Zoom Coupler

907410 16 – 34mm

Light Post Adaptor Sleeve

906135 Storz/Olympus
906137 ACMI
906141 Wolf/Dyonics

Lightsource Cable Adaptor

906133 Male ACMI

PowerTek™ II Plus

906064 Console, Handpiece, Foot Control

PowerTek™ II Plus Shaver Handpiece

906063

PowerTek™ II Plus Shaver Foot Control

906061

PowerTek™ II Plus Console

906062

PowerPump® LT

906050 Console, Foot Control, Remote

PowerPump® LT Stand

906047

Cart with Mini Boom

906082

Saline Bag Level Sensor Module

906022

This material is intended for the Biomet Sports Medicine Sales Force and surgeons only. It is not intended to be redistributed without the express written consent of Biomet Sports Medicine.

PowerTek™ and PowerPump® are trademarks of Biomet Sports Medicine, Inc.



P.O. Box 587, Warsaw, IN 46581-0587 • 800.348.9500 ext. 1501 • ©2006 Biomet Sports Medicine, Inc. All Rights Reserved
Form No. Y-BMT-1012/123106/K



www.biomet sports medicine.com